

Document Title: ON TRACK Reprocessing Inservice / Competency for GI Endoscopes (EVIS, EXERA, EXERA II)

Precleaning	Demonstrated (✓)	Comments
<ul style="list-style-type: none"> ● Wear appropriate personal protective equipment ● Wipe down the insertion tube with a detergent soaked lint free cloth ● Immerse distal tip in detergent and depress suction valve to aspirate detergent for 30 seconds ● Remove distal tip from detergent solution and depress suction valve to aspirate air for 10 seconds ● Turn OFF suction pump and light source ● Attach air/water channel cleaning adapter and set the light source airflow to HIGH ● Immerse the distal tip in clean water ● Depress air/water channel cleaning adapter and feed water for 30 seconds ● Release the air/water channel cleaning adapter for 10 seconds or more to let air through the channel ● Turn OFF light source ● For CF, GIF and other model endoscopes with auxiliary water feeding channels <ul style="list-style-type: none"> ○ Connect auxiliary water tube (MAJ-855) to auxiliary water inlet ○ Immerse distal end in a container of clean water ○ Use a 30 cc syringe to flush detergent solution through auxiliary water tube until no bubbles exit the distal end ○ Use a 30 cc syringe to flush clean water through auxiliary water channel several times ○ Use a 30 cc syringe to flush air through auxiliary water channel until steady stream of air bubbles exits distal tip ● For JF/TJF model duodenoscopes and other model endoscopes with elevator wire channel <ul style="list-style-type: none"> ○ Attach washing tube (MH-974) to the elevator channel plug ○ Use 5 cc syringe to flush detergent solution through the elevator wire until no bubbles exit the distal tip ○ Use 5 cc syringe to flush water through elevator wire channel several times ○ Use 5 cc syringe to flush air through elevator wire channel until a steady stream of bubbles exits distal tip ● Disconnect all removable and reusable parts from the endoscope ● Confirm that water resistant cap is dry and free of debris and attach water resistant cap ● Transport to reprocessing area in covered container 		

Leakage Testing	Demonstrated	Comments
<ul style="list-style-type: none"> ● Fill a basin with clean water ● Confirm that there is no water inside the leakage tester's connector cap ● Connect leakage tester to air source and confirm that air is being emitted ● Connect leakage tester to the endoscope and confirm bending section inflation ● Completely immerse the endoscope in water ● Observe for 30 seconds while angulating the bending section ● Remove endoscope from the water and turn OFF the air source ● Disconnect leakage tester from the air source ● Wait 30 seconds, or until bending section contracts to its pre-expansion size ● Disconnect leakage tester connector cap from venting connector 		

Document Title: ON TRACK Reprocessing Inservice / Competency for GI Endoscopes (EVIS, EXERA, EXERA II)

Manual Cleaning (cont.)	Demonstrated	Comments
<ul style="list-style-type: none"> ● Remove the endoscope and equipment from the detergent solution ● Immerse endoscope and equipment in clean water and gently agitate to rinse ● Reconnect the connector plug and injection tube ● Use a 30 cc syringe to inject 90 cc of water through each side of the injection tube ● For CF, GIF and other model endoscopes with auxiliary water feeding channels <ul style="list-style-type: none"> ○ Attach a 30 cc syringe to the auxiliary water tube (MAJ-855) and inject 90 cc of clean water ● For JF/TJF model duodenoscopes and other model endoscopes with elevator wire channel: <ul style="list-style-type: none"> ○ Attach a 5cc syringe to the washing tube (MH-974) and flush the elevator wire channel with 5cc of clean water ● For CF, GIF and other model endoscopes with auxiliary water feeding channels <ul style="list-style-type: none"> ○ Use a 30cc syringe to inject 90cc of air through the auxiliary water tube ● For JF/TJF model duodenoscopes and other model endoscopes with elevator wire channel: <ul style="list-style-type: none"> ○ Use a 5cc syringe to flush 10cc of air through the elevator wire channel ● Disconnect cleaning equipment. Reprocess as described in Olympus instructions, "Cleaning, Disinfection and Sterilization Procedures for Removable Parts and Cleaning/Reprocessing Equipment" ● Use a lint-free cloth to dry all external surfaces, channel plug and injection tube 		
	Facilities Using the EFP250 Endoscope Flushing Pump to Replace Manual Flushing <ul style="list-style-type: none"> ● Fill two-liter container with detergent solution prepared as recommended by the manufacturer ● Connect inlet hose to EFP250 ● Connect main channel assembly (BLUE) to EFP 250 ● Press the blue START / STOP button ● For scopes with elevator wire or auxiliary water channel <ul style="list-style-type: none"> ○ Connect Special Channel Assembly (PURPLE) ○ Press the PURPLE START / STOP button ● Discard used detergent ● Rinse and fill containers and endoscope basin with clean water ● Press BLUE and PURPLE (if applicable) START / STOP button(s) ● Place screened end of inlet hose in empty container ● Press BLUE and PURPLE (if applicable) START / STOP button(s) ● Use a clean, lint-free cloth to thoroughly wipe and dry the external surfaces of the endoscope and reprocessing equipment and continue reprocessing 	

Document Title: ON TRACK Reprocessing Inservice / Competency for GI Endoscopes (EVIS, EXERA, EXERA II)

Automated Endoscope Reprocessor (AER) High-Level Disinfection*		Comments
AER type: _____ High-level disinfectant type: _____		
<ul style="list-style-type: none"> • Test the potency of the disinfectant solution according to the manufacturer’s instructions • Inspect the connections according to the AER manufacturer’s instructions • Verify that the proper connector is being used for the endoscope being reprocessed • Attach the endoscope connectors/adapters to the AER and endoscope as per the AER manufacturer’s instructions • Operate the AER according to the AER manufacturer’s instructions • Ensure the endoscope is soaked in disinfectant solution according to the liquid chemical germicide manufacturer’s recommendations for time and temperature • Remove the endoscope promptly after the AER cycle is completed • Perform the terminal steps that the AER does not perform (for ex: alcohol and air purge) 		
*FOR FACILITY INTERNAL USE ONLY! Olympus personnel are unable to demonstrate use of individual manufacturer’s AER		

Manual High-Level Disinfection	Demonstrated	Comments
<ul style="list-style-type: none"> • Fill a basin with disinfectant solution • Test the potency of the disinfectant solution according to the manufacturer’s instructions • Attach the channel plug and injection tube to the endoscope • Completely immerse the endoscope and equipment in disinfectant solution • Use a 30cc syringe to inject 90cc of disinfectant into each side of the injection tube. Confirm that no bubbles exit the distal tip • For CF, GIF and other model endoscopes with auxiliary water feeding channels <ul style="list-style-type: none"> ○ Attach the auxiliary water tube (MAJ-855) ○ Use a 30cc syringe to inject 90cc of disinfectant through the auxiliary water tube • For JF model duodenoscopes and other model endoscopes with elevator wire <ul style="list-style-type: none"> ○ Attach the washing tube (MH-974) ○ Use a 5 cc syringe to flush 10cc of disinfectant solution through the elevator wire channel via the washing tube. Confirm that no bubble exit the distal tip ○ While immersed in disinfectant solution, raise the forceps elevator and flush with disinfectant using a 5 cc syringe • Disconnect all equipment from the endoscope • Remove any bubbles that adhere to the surfaces with a clean lint-free cloth • Soak endoscope and equipment for the time and at the temperature recommended by the disinfectant manufacturer • Reconnect the channel plug and injection tube to the endoscope • Attach 30cc syringe to each port on the injection tube and inject 90cc of air • For CF, GIF and other model endoscopes with auxiliary water feeding channels <ul style="list-style-type: none"> ○ Reattach the auxiliary water tube (MAJ-855) ○ Use 30cc syringe to inject 90cc of air through the auxiliary water tube • For JF/TJF model duodenoscopes and other model endoscopes with elevator wire: <ul style="list-style-type: none"> ○ Reattach the washing tube (MH-974) ○ Use 5cc syringe to flush 10cc of air through the elevator wire channel • Remove the endoscope and equipment from disinfectant • Disconnect equipment from endoscope 		

Document Title: ON TRACK Reprocessing Inservice / Competency for GI Endoscopes (EVIS, EXERA, EXERA II)

<i>Rinsing</i>	<i>Demonstrated</i>	<i>Comments</i>
<ul style="list-style-type: none"> ● Fill a basin with sterile water, filtered water, or tap water ● Completely immerse the endoscope and equipment in the water ● Wipe all external surfaces with a lint-free cloth ● Attach the channel plug and injection tube to the endoscope ● Use a 30cc syringe to inject 90cc of water through each side of the injection tube ● For CF, GIF and other model endoscopes with auxiliary water feeding channels <ul style="list-style-type: none"> ○ Attach the auxiliary water tube (MAJ-855) ○ Use 30cc syringe to inject 90cc of water ● For JF/TJF model duodenoscopes and other model endoscopes with elevator wire channel <ul style="list-style-type: none"> ○ Attach the washing tube (MH-974) ○ Use 5cc syringe to flush 15cc of water through the elevator wire channel ● Remove the endoscope and equipment from the water ● Cover the distal tip with a lint-free cloth ● Use a 30cc syringe to inject 90cc of air through each side of the injection tube ● For CF, GIF and other model endoscopes with auxiliary water feeding channels <ul style="list-style-type: none"> ○ Attach the auxiliary water tube (MAJ-855) ○ Use 30cc syringe to inject 90cc of air ● For JF/TJF model duodenoscopes and other model endoscopes with elevator wire channel <ul style="list-style-type: none"> ○ Attach the washing tube (MH-974) ○ Use 5cc syringe to flush 15cc of air through the elevator wire channel ● Disconnect all equipment from endoscope ● Wipe all external surfaces with a lint-free cloth 		

<i>Alcohol Flush</i>	<i>Demonstrated</i>	<i>Comments</i>
<ul style="list-style-type: none"> ● Attach a reprocessed channel plug and injection tube to the endoscope ● Immerse suction port of injection tube in 70% isopropyl or ethyl alcohol ● Attach a 30cc syringe to each port of the injection tube and inject 90cc of alcohol ● Remove the suction port from the alcohol ● For CF, GIF and other model endoscopes with auxiliary water feeding channels <ul style="list-style-type: none"> ○ Attach a reprocessed auxiliary water tube (MAJ-855) ○ Use a 30cc syringe to inject 90cc of alcohol ● For JF model duodenoscopes and other model endoscopes with elevator wire <ul style="list-style-type: none"> ○ Attach a reprocessed washing tube (MH-974) ○ Use 5cc syringe to flush 10cc of alcohol through the elevator wire channel ● Use a 30cc syringe to inject 90cc of air through each side of the injection tube ● For CF, GIF and other model endoscopes with auxiliary water feeding channels <ul style="list-style-type: none"> ○ Attach a reprocessed auxiliary water tube (MAJ-855) ○ Use a 30cc syringe to inject 90cc of air ● For JF model duodenoscopes and other model endoscopes with elevator wire channel <ul style="list-style-type: none"> ○ Attach a reprocessed washing tube (MH-974) ○ Use 5cc syringe to flush 15 cc of air through the elevator wire channel ● Disconnect all equipment from the endoscope ● Use sterile cotton swabs to dry the inside of the air/water and suction cylinders and instrument channel port 		

Document Title: ON TRACK Reprocessing Inservice / Competency for GI Endoscopes (EVIS, EXERA, EXERA II)

<i>Endoscope Storage</i>	<i>Demonstrated</i>	<i>Comments</i>
<ul style="list-style-type: none"> Remove all valves and removable parts from endoscope 		
<ul style="list-style-type: none"> Ensure that angulation locks are in the free position 		
<ul style="list-style-type: none"> For endoscopes with flexible adjustment mechanism, set the insertion tube to maximum flexibility 		
<ul style="list-style-type: none"> Store endoscope in a well-ventilated cabinet 		
<ul style="list-style-type: none"> Hang the endoscope so that the universal cord and insertion tube are hanging vertically and the distal tip of insertion tube is hanging freely 		